

I-SERIES

compatible to: Biomet 3i® / Osseotite® External Hex®



PLATFORM

3,4 mm

4,1 mm

5,0 mm

TITANIUM BASE

for individual milled
Zirconium Abutment
incl. screw
Titanium Grade 5 ELI



I 800

I 810

I 820



SCAN BODY 3D-GUIDE®

for intra + extra
oral scanning



I 9.S3D3.400

I 9.S3D4.150

I 9.S3D4.150



LAB ANALOG

Stainless Steel



I 50

I 51

I 52



STRAIGHT ABUTMENT

Gingiva Height 2,5 mm
incl. screw
Titanium Grade 5 ELI



—

I 110

—



IMPLANT PICK-UP

for open-tray with
integrated screw
Coming soon



I TR-00013.4

I TR-00024.1

—

ABUTMENT SCREW

Recommended tightening
torques 35 Ncm
Titanium Grade 5 ELI



I 61

I 61

I 61



Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

I-SERIES

*compatible to Biomet 3i® /
Osseotite® External Hex®*



4,1 mm



I 10.Mi4.100



I 10.Mi4.100.05



I 5.La4.100



W 12.200.000



I 9.S3D4.150



I 6.Sc4.100



I 61



I-SERIES DIM ANALOG

compatible to: Biomet 3i® / Osseotite® External Hex®



PLATFORM

3,4 mm

4,1 mm

5,0 mm

DIM ANALOG



I 5.DIM.340

I 5.DIM.410

I 5.DIM.500



SLEEVE

for plaster model
Stainless Steel



DIM 5.XSL.A10

DIM 5.XSL.A20

DIM 5.XSL.A30



Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA 